




Ingenio[®] Electroporation Kits & Solution

Enhanced nucleic acid delivery using conventional electroporation devices

- High efficiency electroporation of hard-to-transfect cell lines, stem cells and primary cells
- Compatible with most conventional electroporation devices including Lonza-Amaxa[®], Bio-Rad[®] or Harvard BTX[®]
- Save money and reduce research costs while maximizing results



Ingenio[®] Solution Provides Comparable Efficiency on Amaxa[®] Nucleofector[®] Device.
For experimental details, please visit:
www.mirusbio.com

Delivery by

Mirus.

www.mirusbio.com

Available in Canada from...

MJS
BioLynx
INC.

1-888-593-5969 • www.biolynx.ca • tech@biolynx.ca




Why use Ingenio® Electroporation Kits and Solution?

Ingenio® is a broad spectrum solution that supports high efficiency electroporation with minimal toxicity and replaces standard electroporation solutions including phosphate buffered saline and serum-free media. Ingenio® Kits are compatible with multiple instruments and facilitate a wide range of applications requiring nucleic acid delivery to cells.


Save on Cost Without Compromising Your Results

| Product | Cost/Electroporation* | Savings |
|--|-----------------------|---------|
| Amaxa® Nucleofactor® Kit V (VCA-1003) | \$15.20 | - |
| Ingenio® Electroporation Kit (MIR 50112) | \$9.16 | 40% |


*Based on U.S. list prices from company websites and protocol recommendations (25 reactions)




Ingenio® Solution Provides Comparable Efficiency on the Amaxa® Nucleofactor® Device. Cells were electroporated in parallel with an EGFP reporter vector and assayed at 24 hours post-electroporation by flow cytometry. Two electroporators were used with different electroporation solutions: the Ingenio® Electroporation Kit was used in the Gene Pulser Xcell™ Eukaryotic System (Bio-Rad) and in the Amaxa® Nucleofactor® II/2b Device (Lonza); the Amaxa® Nucleofactor® Kit V was used in the Amaxa® Nucleofactor® II/2b Device, all according to manufacturers' recommendations.




Ingenio® Outperforms Other Electroporation Solutions in Efficiency and Viability. Cells were electroporated in parallel with an EGFP reporter vector using either Ingenio® Electroporation Solution, PBS or Gene Pulser® Electroporation Buffer (Bio-Rad) on the Gene Pulser Xcell™ Eukaryotic System. EGFP expressing cells were identified 24-hours post-electroporation by flow cytometry and presented as a percentage of live cell population.








Efficient Plasmid DNA Delivery in Many Cell Types Using the Amaxa® Nucleofector® Device. Cells were assayed at 24 hours post-electroporation by flow cytometry and reported as percentage of live cell population. Visit www.mirusbio.com for ideal pulse conditions. (*Primary cell types)



Efficient Plasmid DNA Delivery in Many Cell Types Using the Bio-Rad® GenePulser Xcell™ System. Cells were assayed at 24 hours post-electroporation by flow cytometry and reported as percentage of live cell population. Visit www.mirusbio.com for ideal pulse conditions. (*Primary cell types)



High Efficiency Plasmid DNA Electroporation of Human Induced Pluripotent Stem (iPS) Cells using Ingenio®. The Ingenio® Electroporation Kit was used to transfect iPS cells on the Amaxa® Nucleofector® II/2b Device. Cells were electroporated with ZsGreen expressing plasmid (Clontech) and visualized 24 hours post-transfection and imaged under 4X objective with an Olympus IX71® Inverted Microscope. Cells were also assayed 24 hours post-transfection on an Accuri® Cytometer.

| PRODUCT | PRODUCT NO. | QUANTITY |
|--|-------------|------------------------|
| Ingenio® Electroporation Kits (solution, 0.2 cm cuvettes, cell droppers)  Compatible with Lonza-Amaxa® Nucleofector® II/2b devices | MIR 50112 | 25 RXN |
| | MIR 50115 | 50 RXN |
| | MIR 50118 | 100 RXN |
| Ingenio® Electroporation Kits (solution, 0.4 cm cuvettes, cell droppers)  Compatible with Bio-Rad® and Harvard-BTX® devices | MIR 50113 | 25 RXN |
| | MIR 50116 | 50 RXN |
| | MIR 50119 | 100 RXN |
| Ingenio® Electroporation Solution  | MIR 50111 | 25 RXN (6.25 ml) |
| | MIR 50114 | 50 RXN (12.5 ml) |
| | MIR 50117 | 100 RXN (25 ml) |
| Ingenio® Electroporation Accessories   | MIR 50120 | 0.2 cm cuvettes (25PK) |
| | MIR 50121 | 0.2 cm cuvettes (50PK) |
| | MIR 50122 | 0.4 cm cuvettes (25PK) |
| | MIR 50123 | 0.4 cm cuvettes (50PK) |
| | MIR 50124 | Cell Droppers (25PK) |
| | MIR 50125 | Cell Droppers (50PK) |

Prove it to Yourself with a **FREE SAMPLE**



Reagent Agent®

Transfection reagent recommendations based on citations, customer feedback, and in-house transfection data. Find the ideal delivery solution for your experiment: www.mirusbio.com/RA

All trademarks are the property of their respective owners. Ingenio and Mirus Bio are registered trademarks and Reagent Agent is a service mark of Mirus Bio LLC. Amaxa Nucleofector II/2b Device is a registered trademark of Lonza Group Ltd. Bio-Rad and Gene Pulser are registered trademarks and Xcell is a trademark of Bio-Rad Laboratories, Inc. BTX is a registered trademark of Harvard Apparatus, Inc. Accuri is a registered trademark of Becton, Dickinson and Company. IX71 is a registered trademark of Olympus.