

Exosomes

Purification | Detection
Tracking | NTA Service

Exosome Experts

Available in Canada from...

MJS
BioLynx
INC.

1-888-593-5969 • www.biolynx.ca • tech@biolynx.ca



Purification


Exo-spin™ Exosome Purification Kit



An independently validated, flexible system for isolation of exosomes from a variety of sources

- Excellent yields and high levels of purity**
 Exosomes with ultra-low protein and rRNA contamination.
- No ultracentrifugation required**
 Protocol provides consistent results every time.
- Simple and reliable**
 Isolate intact whole exosomes for functional studies.

Exo-spin™ Exosome Isolation Protocol



Step 1: Remove cells and cellular debris.


Step 2: Use Exo-spin™ Exosome Precipitation Buffer to precipitate exosomes.

Step 3: Add the exosome-containing pellet to the Exo-spin™ columns for Size Exclusion Chromatography (SEC) purification.

| | Exo-spin™ | Precipitation | Ultracentrifugation |
|---------------|-----------|---------------|---------------------|
| Yield | High | Medium | Medium |
| Purity | High | Low | Medium |
| Speed | Quick | Quick | Slow |
| Cost | Low | Medium | Medium |

Exosome Extraction Profiles

Exo-spin™ Midi Columns allow for gravity assisted fractionation, where exosomes are separated from the vast majority of proteins.



Exosomes were isolated using Exo-spin™ Midi Columns (Cat Code EX04) from 60 ml of conditioned medium generated by a human breast carcinoma cell line. Fractions (each of 500 µl) were collected and analysed to (1) evaluate particle numbers and (2) measure absorbance at 280 nm to evaluate protein concentration.

| Cat code | Product name |
|----------|------------------------|
| EX01 | Exo-spin™ |
| EX02 | Exo-spin™ blood |
| EX03 | Exo-spin™ mini columns |
| EX04 | Exo-spin™ midi columns |

Detection


TRIFic™ Exosome Detection Kit



An exquisitely sensitive Europium Time-Resolved Immunofluorescence assay for exosome markers

- Simplicity**
 Assay is clear and simple allowing for high reproducibility.
- Sensitivity**
 Europium-labeled antibodies are used for detection enabling the use of TRF.
- Specificity**
 Only antigens displayed in multiple copies are detected.

TRIFic™ Exosome Assay




Step 1: Biotinylated antibody is bound to streptavidin coated assay plates.


Step 2: Biological samples are added. Exosomes and any free antigen are captured by the antibody.

Step 3: Europium-labeled antibody of the same clone as the one used in Step 1 is added and binds specifically to exosome antigen. The epitopes of bound monomers are already occupied and not detected. Samples are read on a time-resolved fluorescence plate reader.

Sample profiling



CD9 TRIFic™ exosome assay performed for 19 different cell lines. Off the scale readings indicated (x1000) on individual bars. Samples generating vastly different signal intensities can be measured due to the broad signal range covered by the assay.



TRIFic™ exosome assay analysis of 20 urine samples shows great variation in CD9 content between samples.

| Cat code | Product name |
|----------|--------------------------------------|
| EX101 | TRIFic™ CD9 Exosome Assay, 96 wells |
| EX102 | TRIFic™ CD63 Exosome Assay, 96 wells |
| EX103 | TRIFic™ CD81 Exosome Assay, 96 wells |
| EX-P31 | Wash Buffer (Concentrate 25x), 20 ml |

Tracking


ExoFLARE™ Exosome Tracking Kit




An exosome tracking system using protein tags that modulate the fluorescence of pro-fluorophores

- High signal intensity**
 The signal is brighter and longer than for standard fluorescent proteins.
- Low background fluorescence**
 The rapid turnover of the dye allows for very low photobleaching.
- Non-cytotoxic**
 Further experiments can be performed with same exosomes.

ExoFLARE™ Exosome Tracking System



In vitro exosome imaging



ExoFLARE™ constructs were transiently transfected into DU145N cells. ExoFLARE™ cell permeable dye was added to the media and cells were imaged using a confocal fluorescence microscope.

Red = ExoFLARE™ tagged protein;
Green = Hoechst.

| Cat code | Product name |
|----------|----------------------------------|
| EX301 | ExoFLARE™ CD9 (impermeable dye) |
| EX302 | ExoFLARE™ CD63 (impermeable dye) |
| EX303 | ExoFLARE™ CD81 (impermeable dye) |
| EX304 | ExoFLARE™ CD9 (permeable dye) |
| EX305 | ExoFLARE™ CD63 (permeable dye) |
| EX306 | ExoFLARE™ CD81 (permeable dye) |
| EX401 | EF Red xc cell impermeable dye |
| EX402 | EF Red s cell permeable dye |

NTA Profiling


Nanoparticle Tracking Analysis service provided using ZetaView® Instrument




Exosome characterization service for analysis of particle size and particle concentration

- High quality**
 The service is performed in our labs by highly qualified scientists.
- Quick turnaround times**
 Full reports are e-mailed within 5 to 10 business days.
- Competitive price**
 Very low price per sample without the need of purchasing the equipment.

Example of an NTA report generated with the Zetaview®.



Zetaview®
Instrument.



Individual particle movement is tracked and recorded for characterization of exosome samples.

Image provided by Particle Matrix GmbH.

| Cat code | Product name |
|----------|--|
| ZV-1 | Analysis service set up (purchase one per order) |
| ZV-12 | NTA analysis of a single sample |

Cell Guidance Systems' reagents and services enable control, manipulation and monitoring of the cell, both *in vitro* and *in vivo*

Growth Factors

- Recombinant
- Sustained Release

Exosomes

- Purification
- Detection
- Tracking
- NTA Service

Small Molecules

Cell Counting Reagent

Matrix Proteins

Cell Culture Media


- Pluripotent Stem Cells
- Photostable
- *In Vitro* Blastocyst Culture
- ETS-embryo Culture
- Custom Manufacturing Service

Gene Knock-Up System

Cytogenetics Analysis



www.cellgs.com



1-888-593-5969 • www.biolynx.ca • tech@biolynx.ca