

Sample Thawing Stations

"Developed by a scientist at a major genomic diagnostic company to improve productivity and quality of analytical results"

New HeliPad accessory racks expand the versatility and productivity of the Heliport Thawing Stations



The new HeliPad racks can be placed anywhere on the Heliport deck and can be mixed with other rack consumables.

The small HeliPad rack has 40 holes for 2 ml size tubes (11mm).

The large HeliPad rack has 14 medium size holes (17mm) for 15ml size tubes and 10 large holes (30mm) for large conicals.

Place anywhere on the deck for maximum versatility.



The Heliport can accommodate up to 5 HeliPad racks. Above configuration: 120 2ml tubes, 30 15 ml tubes and 20 large conicals

Ordering Info: BOX 130080000 - Helipad Rack Set BOX 130081000 - Helipad Small Tube BOX 130082000 - Helipad Large Tube BOX 130070000 - Heliport Thawing BOX 130071000 - Heated Heliport Plus



1-888-593-5969 • www.biolynx.ca tech@chromspec.com





The HeliPads simply snap into

place anywhere on the

Heliport deck.