


USING THE CRYOTRAY

The Seal'N Freeze Cryotray was developed for use in our lab to address common challenges in frozen tissue sample processing and storage.



The Cryotray can be used with the O.C.T. compound of your choosing.

To close the mold for freezing, storage, and transport, **press firmly on opposite sides of the lid until it snaps into place.**

The Cryotray is designed for **easy separation of each mold** - carefully fold and tear away from the center connector to separate each sample.

Each mold is constructed of the highest-quality, most durable plastic for **smooth and easy extraction of each frozen block** for cryosectioning.



VISIT www.sealnfreeze.com

**SEAL'N
FREEZE**
A VANESSA RESEARCH BRAND



CRYOTRAY[®]
QUALITY. PRECISION. EFFICIENCY

The highest standard in frozen tissue sample processing and storage

- ✓ Sealable lids for optimal protection
- ✓ Perfect, uniform blocks every time
- ✓ Higher quality, lower costs
- ✓ Freeze four samples simultaneously

Available in Canada from...

MJS
BioLynx
INC.

1-888-593-5969 • www.biolynx.ca • tech@biolynx.ca



FOR THE SCIENTIST



OPTIMAL SAMPLE PROTECTION

Each mold includes a **sealable lid** for optimal sample protection and storage.



PERFECT FROZEN SAMPLES

Trays are constructed of high-quality, durable plastic for **smooth, easy extraction** of your frozen sample from the mold.



BATCH PROCESSING

The four mold Cryotray allows for batch processing of tissue samples. Level Cryotray quadrants create the ideal environment for **precision and accuracy**.



CONVENIENT + HANDS-FREE

Place the Cryotray in the Seal'N Freeze box, close the lid, and wait - no need to manually hold each sample above your freezing media.



The Seal'N Freeze Cryotray features conveniently attached lids

FOR THE LAB MANAGER

THE MOST COST-EFFECTIVE MOLD

For just \$0.37 per mold, Seal'N Freeze offers significant time-savings on sample freezing, along with a **free cryopreservation box** for hands-free processing.

75% INCREASE IN EFFICIENCY

FREEZE TIME PER SAMPLE, USING LIQUID NITROGEN



WITHOUT SEAL'N FREEZE



WITH SEAL'N FREEZE

**Per sample freeze-time calculated based on time to freeze one tray, divided by 4*

Decrease the time your lab spends on freezing by 75% - allowing your team to spend time on more important work.

PRODUCTS

THE SEAL'N FREEZE CRYOTRAY

50 Pack of
Cryotrays
(200 molds)



Mold dimensions:
15L x 15W x 5D mm (4 molds per tray)

THE SEAL'N FREEZE BOX



Maximize your efficiency by using the Cryotray with our compact Seal'N Freeze Box. Simply fill the box with your choice of freezing media and place the Cryotray on the platform for convenient, hands-free freezing.

MJS
BioLynx
INC.

1-888-593-5969

www.biolynx.ca • tech@biolynx.ca