



bringing innate immunity to the **next level** 

Available in Canada from...



1-888-593-5969 • www.biolynx.ca • tech@biolynx.ca

### Oxidative stress

Hycult Biotech is a leading world-class manufacturer of research reagents in the field of innate immunity. We are a specialized partner in the development and manufacturing of assays, antibodies, proteins and peptides. With the focus on innate immunity, we pay special attention to products relating to antimicrobial peptides, complement, acute phase proteins, TLR, LPS and microbial toxins.

Oxidative stress is involved in many human diseases, such as atherosclerosis, Parkinson's disease and Alzheimer's disease. In addition, it is assumed to be important in ageing. However, reactive oxygen species can also be beneficial, as they are used by the immune system to attack and kill pathogens. Detection or measurement of oxidation markers is helpful to assess oxidant activity and to monitor the effectiveness of the antioxidant system. Superoxides and other free radicals contribute, together with many other factors like inflammation, radiation, and carcinogen exposure, to DNA damage. Antibodies and assays for detection of DNA adducts, such as ethenoadenosine, BPDE-DNA, and 8-oxoguanine, are valuable markers for studies on DNA damage. Nitrotyrosine has been identified as a marker of inflammation and production of NO. Various pathways including the formation of peroxynitrite lead to nitrotyrosine production. Since nitrotyrosine is a stable end product of peroxynitrite oxidation, assessment of its plasma concentration is a useful marker of NO-dependent damage in vivo, taking into account that most proteins have a longer half life in the circulation than nitrite/nitrate (NOX). Nitrotyrosine has been detected by immunohistology in various inflammatory processes, including atherosclerotic plaques, celiac disease, rheumatoid arthritis, chronic renal failure and septic shock.

#### Validated ELISA kit sample types:

BF:	Bronchoalveolar fluid	CS:	Cell culture supernatant
F:	Feces	O:	Other biological fluids
P:	Plasma	S:	Serum
SP:	Sputum	TH:	Tissue Homogenate
U:	Urine		

ELISA kits				
Product name	Species	Sample types	Cat. #	
Arginase I	Human	CSPS	HK322	
MGO, Competitive	N/A	FP	HIT503	
MPO	Human	CS P SP	HK324	
MPO	Mouse	CS P TH	HK210	
MPO	Rat	CS P TH	HK105	
Nitrotyrosine	N/A	PO	HK501	
Nitrotyrosine, Competitive	N/A	FPS	HIT501	


Proteins				
Product name	Quantity	Cat. #		
Oxidized PAPC	1 mg	HC4035		
Oxidized PAPC	5 mg	HC4036		
PAPC	5 mg	HC4043		

## Can't find the sample type you want for an ELISA kit?

Ask the support team to assist you with your research needs on <a href="mailto:support@hycultbiotech.com">support@hycultbiotech.com</a>!

## How many tests would you like?

Almost all our ELISA kits are available in both a 1-plate and a 2-plate format!




Nitrotyrosine in human liver. Staining of paraffin section with clone HM.11 (Cat. # HM5001)




### Validated antibody applications:

F:	Frozen sections	FC:	Flow cytometry
FS:	Functional studies	IA:	Immunoassay
IF:	Immunofluorescence	IP:	Immunoprecipitation
P:	Paraffin sections	W:	Western blot

Antibodies					
Product name	Clone	Species	Applications	Cat. #	
Arginase I	6G3	Human	F FC IA	HM2162	
Arginase I	9C5	Human	F IA IP	HM2163	
Arginase I	D1.2	Human	FC IA W	HM2324	
BPDE	8E11	N/A	IA	HM5008	
BPDE-DNA	5D11	N/A	FC IA IF IP P	HM5007	
Caveolin-1	7C8	Rat	IF IP W	HM3014	
Chlorotyrosine	pAb	N/A	FIAP	HP5002	
CML	CML26	Human	FIAIFPW	HM5013	
Ethenoadenosine	1G4	N/A	F FC IA IP W	HM5005	
Halotyrosine	BTK-94C	N/A	IA IF IP W	HM5016	
LOX-1	23C11	Human	FC FS IF P	HM2138	
MG-H1	1H7G5	N/A	IA IF W	HM5017	
MGO-modified proteins	MGO-1	Human	IAPW	HM5014	
MPO	266-6K1	Human	IA W	HM2164	
MPO	B7	Human	F FC IA W	HM2288	
MPO	pAb	Human	FIAP	HP9048	
MPO	8F4	Mouse	F FC IA IF	HM1051	
MPO	B7	Mouse	F FC IA W	HM1135	
MPO	2D4	Rat	FC IA	HM3030	
MPO	B7	Rat	F FC IA W	HM3039	
Nitrotyrosine	HM.11	N/A	FIAPW	HM5001	



LOX-1 in human spleen. Staining of paraffin section with clone 23C11 (Cat. # HM2138)



Arginase I in human liver. Staining of frozen section with clone 6G3 (Cat. # HM2162)  $\,$ 





# What we stand for

· Original manufacturer

Innate immunity focus

• > 20 years experience

- Reliable and robust
- Support within one working day
- Scientific drive

Hycult Biotech, bringing innate immunity to the next level



www.hycultbiotech.com

Hycult Biotech is distributed by:

Even though the information provided is verified and believed to be accurate, there is always a chance that it is imprecise or may be outdated. Please make sure to visit our website for the most recent product information or contact us in case you have any questions or comments.



1-888-593-5969 www.biolynx.ca • tech@biolynx.ca