



Station One

Thaw Right!

The first of its kind! The Station One thaws frozen reagents and samples in a fraction of the time using only ambient air. Its fixed platform enables high throughput sample thawing that is repeatable, proceduralizable, and safe for even the most delicate samples and reagents.

The Station One channels ambient air over and around samples. The result is up to a 75% decrease in sample and reagent thawing times without any application of heat. Included accessory racks enable thawing of most industry standard tubes from 1 to 50ml.

The Station One's large thawing surface can accommodate two microplates or 2D barcoded tube racks at a time. There's ample surface area for thawing media bottles too.

The Station One unit has a compact footprint and consumes minimal power.



What is the cost of your time?

Our thaw stations can save you time and money while assuring you consistent inputs into your manual and automated laboratory processes. Microplates, tubes and even 2D barcoded tube racks can be ready for use in a fraction of the time. Our surface capacity and interchangeable racks provide you maximum speed and flexibility. There is no single product on the market that can offer all of this.



Achieve Equilibrium with the BOX Station One!

Product Number: BOX13000725

Box Scientific
www.boxscientific.com

Available in Canada from...



1-888-593-5969 • www.biolyx.ca • tech@biolyx.ca

